

DATE \_\_\_\_\_ REQUESTED RETURN DATE \_\_\_\_\_

DR. \_\_\_\_\_

AGE \_\_\_\_\_  
 MALE  
 FEMALE  
 PLEASE CALL

PATIENT'S NAME \_\_\_\_\_

PLEASE RETURN CASE  
 BISQUE BAKE     GLAZED     FOR TRY-IN  
 CASTING     ASSEMBLED FRAMEWORK

**P F M**

**ALLOY**  
 HIGH-NOBLE (precious)     WHITE or  YELLOW  
 NOBLE (semi-precious)    (CIRCLE ONE)  
 CAPTEK     BASE ALLOY (non-precious)

**If insufficient room:**  
 REDUCE AND MARK     PLEASE CALL     REDUCTION COPING

**PONTIC DESIGN**  
 SANITARY     FULL RIDGE     MODIFIED     BULLET

PORCELAIN BUTT JOINT:  
 180°  
 360°

**MARGIN DESIGN**  
 BAND AT BUCCAL:     Hair Line     1 mm     2 mm     No metal exposed  
 BAND AT LINGUAL:     Hair Line     1 mm     2 mm     No metal exposed

FABRICATION PURPOSE:     Correct Malalignment     Close Spaces  
     Increase Length     Color Change

SURFACE TEXTURE:     Smooth     Moderate     Heavy

Rx \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

A 1.5% PER MONTH SERVICE CHARGE WILL BE APPLIED TO ALL PAST DUE BALANCES. EACH PRESCRIPTION MUST BE COMPLETED & SIGNED


DR. SIGNATURE \_\_\_\_\_ LIC. NO. \_\_\_\_\_

PLEASE SEND MORE RX'S  
 PLEASE SEND MORE BAGS

*Thank you!*

**SHADING CHART**

Shade of prepared teeth: \_\_\_\_\_  
 Shade desired: \_\_\_\_\_



**GOLDEN PROPORTION:**  
 Central Incisor Length \_\_\_\_\_ mm  
 Width \_\_\_\_\_ mm

**VALUE :**  
 High (bright)  
 Medium  
 Low

**OCCLUSAL STAIN:**  
 None  
 Light  
 Medium  
 Dark

**METAL FREE ■ COSMETIC**

**LAMINATES**  
 CERAMCO     NORITAKE

**ALL CERAMIC CROWNS & BRIDGES**  
 E-MAX     ZIRCONIA CROWN  
 BRUXZIR     ZIRCONIA BRIDGE  
 ALL PORCELAIN CROWN

**INLAY-ONLAY**  
 IPS EMPRESS MILLED CERAMIC  
 E-MAX PRESSABLE CERAMIC  
 PREMISE INDIRECT COMPOSITE

**DESIGN YOUR CASE HERE:**  
 STUDY MODELS SHOULD COME WITH ALL CASES

